



## Technical details

---



No. of outputs:	1
OUTPUT VOLTAGE:	56V
WAVE SHAPE:	SINUSOIDAL
WIDTH OF IMPULSES:	140 mSec
NO. OF PROGRAMS:	1
OUTPUT CURRENT :	56 mA (load 1kohm)MAX
CARRIER FREQUENCY:	7.500 Hz
MODULATION FREQUENCY:	70 Hz
SUPPLY VOLTAGE:	9V (ALKALINE BATTERY)

CONTRACTION TIME: 3 Sec.

RECOVERY TIME: 3 Sec.

MARKING:



## Tua Viso is safe

Tua viso is safe and complies with the most stringent international regulations. It is entirely designed and manufactured by Vupiesse Italia, a company with a Certified Quality System in conformity with **ISO 9001:2000 and ISO13485:2004**, issued by the German institutes **TÜV RHEINLAND**. The company has considerable experience in the manufacture of professional electro-medical equipment, used mainly by the major hospitals and centres that deal with sports medicine, physiotherapy and beauty medicine.



---

## Precautions and contraindications

- During the menstrual period, skin conduction changes; stimulation could prove less comfortable.
- Stimulation causes reddening of the skin where the electrodes are placed, especially in the case of persons with sensitive skin. Such reddening is perfectly normal and usually disappears quickly.
- Do not apply the electrodes to nevi for prolonged time.
- Do not use on mouth and eyes.
- Do not use on the abdomen of expectant mothers.
- Do not use Tua Viso on children, people with pace-makers or those in poor psychophysical conditions or with physiological or organic disorders.

IN CASE OF DOUBT, SEE A DOCTOR.